






Afname- en verzendmedia: Micro afnamebuisen

REF 365979 Microgelbuis $\geq 400 - \leq 600 \mu\text{l}$ Microtainer	Geel	
REF 365975 K₂-EDTA $\geq 250 - \leq 500 \mu\text{l}$ Microtainer	Lila	
REF 366988 Lithium heparine $\geq 400 - \leq 600 \mu\text{l}$ Microtainer	Lichtgroen	
REF 365966 Lithium heparine $\geq 50 - \leq 500 \mu\text{l}$ Microtainer	(donker)Groen	
REF 41.1506.005 Citraatbuis 1300 μl 130 μl citraatsolution Sarstedt	(licht)Groen	
REF 365993 NaF microbuis $\geq 400 - \leq 600 \mu\text{l}$ Microtainer	Grijs	
REF 365964 Stolbuis (zonder toevoegingen) o.a. voor medicijnspiegel Microtainer	Rood	

Voor meer informatie over aanvragen, zie "Laboratorium Informatie Gids": via www.labgids.nl.

Voor meer informatie over bloedafname d.m.v. micro afname zie: /BLO/WI103019 "[Bloedafname d.m.v. micro-afname](#)".

Voor meer informatie over minimale hoeveelheden zie: /BLO/BL103013 "[Bloedafname minimale hoeveelheden](#)".

Versie: 8, augustus 2021